

Product Datasheet

ADAR1 Rabbit mAb Antibody (orb1950638)

Description	ADAR1 Rabbit mAb Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC, WB
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P55265
Expiration Date	12 months from date of receipt.