

Product Datasheet

PCAF Rabbit mAb Antibody (orb1950727)

Description	PCAF Rabbit mAb Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IP, WB
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q92831
Expiration Date	12 months from date of receipt.