

Product Datasheet

HXB2 antibody (orb19513)

Description

HXB2 antibody

Species/Host

Rabbit

Reactivity

Human, Rat

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications

Immunogen

KLH conjugated synthetic peptide derived from human HXB2. Please contact us for the exact immunogen sequence. The peptide is available as orb393972.

Target

HXB2

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW

38 kDa

Uniprot ID

P10683

NCBI

015973, 057057

Entrez


51083

Dilution Range


IHC-P: 1:100-800, WB: 1:200-800

Expiration Date


12 months from date of receipt.



IHC-P image of human lung cancer tissue ...



Immunohistochemical staining of rat spin...



Immunohistochemical staining of human lu...