

Product Datasheet

PAX6 Rabbit mAb Antibody (orb1951361)

Description	PAX6 Rabbit mAb Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IHC, WB
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P26367
Expiration Date	12 months from date of receipt.