

Product Datasheet

STAG2 antibody (orb19515)

Description	Goat polyclonal to STAG2
Species/Host	Goat
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	STAG2
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	This antibody is expected to recognise isoform a (NP_001036214.1 and NP_001036215.1) and isoform b (NP_001036216.1 and NP_006594.3).ELISA: Peptide ELISA: antibody detection limit dilution 1:128000.WB: Approx 150kDa band observed in lysates of cell line HeLa (calculated MW of 146kDa according to NP_001036214.1 and NP_001036215.1). Recommended concentration: 0.01 µg/ml.
Protein Sequence	REPKRLRPEDS
Clonality	Polyclonal
MW	146; 146; 141; 141
Entrez	10735
Dilution Range	ELISA: 1:128000, WB: 0.01 µg/ml
Expiration Date	12 months from date of receipt.



orb19515
(0.01µg/ml)
staining of
HeLa ce...