

Product Datasheet

PLC gamma 1 Rabbit mAb Antibody (orb1951581)

Description	PLC gamma 1 Rabbit mAb Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IP, WB
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P19174
Expiration Date	12 months from date of receipt.