

Product Datasheet

Histone H3 (acetyl K14) Rabbit mAb Antibody (orb1951813)

Description	Histone H3 (acetyl K14) Rabbit mAb Antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IHC, IP, WB
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P68431
Expiration Date	12 months from date of receipt.