

## Product Datasheet

### alpha Tubulin Rabbit mAb Antibody (orb1951998)

|                            |                                                                                |
|----------------------------|--------------------------------------------------------------------------------|
| <b>Description</b>         | alpha Tubulin Rabbit mAb Antibody                                              |
| <b>Species/Host</b>        | Rabbit                                                                         |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                              |
| <b>Conjugation</b>         | Unconjugated                                                                   |
| <b>Tested Applications</b> | FC, ICC, IHC, IP, WB                                                           |
| <b>Storage</b>             | Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle. |
| <b>Note</b>                | For research use only                                                          |
| <b>Clonality</b>           | Recombinant                                                                    |
| <b>Uniprot ID</b>          | <b>P68366</b>                                                                  |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                |