

Product Datasheet

alpha Tubulin Rabbit mAb Antibody (orb1951999)

Description	alpha Tubulin Rabbit mAb Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC, IP, WB
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P68366
Expiration Date	12 months from date of receipt.