

Product Datasheet

Cyclooxygenase 2 Rabbit mAb Antibody (orb1952166)

Description	Cyclooxygenase 2 Rabbit mAb Antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IP, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q05769
Expiration Date	12 months from date of receipt.