

Product Datasheet

cAMP Protein Kinase Catalytic Subunit Rabbit mAb Antibody (orb1952179)

Description	cAMP Protein Kinase Catalytic Subunit Rabbit mAb Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IHC, IP, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P17612
Expiration Date	12 months from date of receipt.