

Product Datasheet

Monkey Immunoglobulin G (IgG) ELISA Kit (orb1952282)

Description	Monkey Immunoglobulin G (IgG) ELISA Kit
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 3.12-200 µg/mL. Sensitivity: 1.36 µg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 50 µL. Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Expiration Date	Please enquire.