

Product Datasheet

Monkey N-terminal pro-brain natriuretic peptide (NT-proBNP) ELISA Kit (orb1952287)

Description Monkey N-terminal pro-brain natriuretic peptide (NT-proBNP) ELISA Kit

Tested Applications ELISA

Note For research use only

Application notes Detect Range: 1.56-100 ng/mL. Sensitivity: 0.62 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 50 µL. Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Expiration Date Please enquire.