

Product Datasheet

Human Growth Differentiation Factor 15 (GDF15) ELISA Kit (orb1952291)

Description	Human Growth Differentiation Factor 15 (GDF15) ELISA Kit
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 78.1-5000 pg/mL. Sensitivity: 31 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 50 µL. Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Expiration Date	Please enquire.