

Product Datasheet

Human Interferon Gamma (IFN γ) ELISA Kit (orb1952292)

Description	Human Interferon Gamma (IFN γ) ELISA Kit
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 31.2-2000 pg/mL. Sensitivity: 11.74 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 50 μ L. Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Expiration Date	Please enquire.