


Product Datasheet

TGM2 antibody (orb19526)

Description	Goat polyclonal antibody to TGM2
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	TGM2
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:128000. WB: Approx 90kDa band observed in Human Lung lysates (calculated MW of 77.3kDa according to NP_004604.2). Recommended concentration: 0.03-0.1 µg/ml.
Protein Sequence	NFESDKLKAVKGFR
Clonality	Polyclonal
MW	77.3
Entrez	7052
Dilution Range	ELISA: 1:128000, WB: 0.03-0.1 µg/ml, IHC-P: 5-10 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Lung lysa...



Western blot analysis of A549 lysate usi...