



Product Datasheet

TCR gamma/delta Antibody (APC) (orb1952611)


Description	The antigen-specific T cell receptor (TCR) is composed of either alpha and beta subunit, or gamma an...
Reactivity	Human, Primate
Conjugation	APC
Tested Applications	FC
Target	TCR gamma/delta
Preservatives	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: The reagent is designed for analysis of human blood cells using 10 µl reagent / 100 µl of whole blood or 10 ⁶ cells in a suspension. The content of a vial (1 ml) is sufficient for 100 tests.
Isotype	Mouse IgG1 kappa
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.



Flow
cytometry
multicolor
surface
staini...



Separation
of human
CD3
positive TCR
gam...



Flow
cytometry
surface
staining
pattern ...