

Product Datasheet

Recombinant SARS-CoV-2 Nucleoprotein (orb1952703)

Description

Recombinant SARS-CoV-2

Reactivity

Virus

Tested

ELISA, WB

Applications
Preservatives

Lyophilized from a 0.2 µm filtered solution in 10 mM Hepes, 500 mM NaCl with 5% trehalose, pH 7.4

Storage

The lyophilized protein is stable at -20°C for up to 1 year, the protein solution can be stored at 2-8°C for up to 1 week. For extended storage, it is recommended to further dilute in a buffer containing 0.1% BSA (make sure BSA will NOT interfere with your experiment) and store in working aliquots at -20°C to -80°C for 6 months. Avoid repeated freeze/thaw cycle.

Tag

N-terminal His-EGFP Tag

Note

For research use only

Application notes

Centrifuge the vial prior to opening, reconstitute in sterile distilled water or buffer of interest to concentration of 0.1-1mg/ml by gently pipetting 2-3 times, don't vortex

Purity

Greater than 95% as determined by SDS-PAGE

Source

E.Coli

MW

76.7 kDa

Uniprot ID
P0DTC9
Expiration Date

6 months from date of receipt.