



Product Datasheet

PPARGC1A antibody (orb19529)


Description	Goat polyclonal antibody to PPARGC1A
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, IHC
Target	PGC1A/PPARGC1A
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:64000.WB: Approx 100kDa band observed in Rat Stomach lysates (calculated MW of 91.0kDa according to NP_037393.1. Recommended concentration: 1-3 µg/ml.
Protein Sequence	DGLFDDSEDESDK
Clonality	Polyclonal
MW	91.0; 91.7; 89.6; 77.1; 91.5
Entrez	10891
Dilution Range	ELISA: 1:32000, WB: 0.5-2 µg/ml, IF/ICC: 10µg/ml, FACS: 10ug/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Rat Stomach lys...



Immunofluorescence analysis of paraforma...



Immunofluorescence analysis of paraforma...