

Product Datasheet

Mouse Anti-Human CD3/CD19/CD45, FITC/PE/PE-Cy5 Conjugated mAb [Out of stock] (orb195290)

Description	Mouse Anti-Human CD3/CD19/CD45, FITC/PE/PE-Cy5 Conjugated mAb
Species/Host	Mouse
Conjugation	PE/Cy5
Tested Applications	FC
Storage	Store at 4°C.DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.
Note	For research use only
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.