

Product Datasheet

Human SPINT1/HAI-1 (Serine Peptidase Inhibitor, Kunitz Type 1) ELISA Kit (orb1952966)

Description	Human SPINT1/HAI-1 (Serine Peptidase Inhibitor, Kunitz Type 1) ELISA Kit
Reactivity	Human
Range	31.2-2000pg/ml
Target	SPINT1/HAI-1
Storage	2-8 °C for 6 months
Note	For research use only
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.