

## Product Datasheet

### Mouse anti-Human CD45, Purified mAb [Out of stock] (orb195300)

|                            |                                                          |
|----------------------------|----------------------------------------------------------|
| <b>Description</b>         | Mouse anti-Human CD45,<br>Purified mAb                   |
| <b>Species/Host</b>        | Mouse                                                    |
| <b>Tested Applications</b> | FC                                                       |
| <b>Storage</b>             | Store at 4°C.DO NOT FREEZE.<br>LIGHT SENSITIVE MATERIAL. |
| <b>Note</b>                | For research use only                                    |
| <b>Clonality</b>           | Monoclonal                                               |
| <b>Expiration Date</b>     | 12 months from date of receipt.                          |