

Product Datasheet

Mouse IgG1 Isotype Control, PE Conjugated mAb [Out of stock] (orb195302)

Description	Mouse IgG1 Isotype Control, PE Conjugated mAb Antibody
Species/Host	Mouse
Conjugation	PE
Tested Applications	FC
Form/Appearance	Phosphate-buffered solution, pH 7.4, containing 0.09% sodium azide and 0.2% (w/v) rAlbumin
Storage	Store at 4°C.DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.
Note	For research use only
Application notes	Vol.per.Test: 10 µl/Test
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.