

Product Datasheet

Robo3 antibody (orb195324)

Description	Goat polyclonal antibody to Robo3
Species/Host	Goat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA
Target	Robo3.1 (mouse aa1338-1351)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Preliminary experiments gave bands at approx 100kDa and 80kDa in Mouse fetal Brain lysates after 1µg/ml antibody staining. Please note that currently we cannot find an explanation in the literature for the bands we observe given the calculated size of 150
Protein Sequence	PYGRPSFLSHGQGT
Clonality	Polyclonal
MW	150
Dilution Range	ELISA: 1:128000, WB: 1 µg/ml
Expiration Date	12 months from date of receipt.