

Product Datasheet

PHF19 antibody (orb195332)

Description	Goat polyclonal antibody to PHF19
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Target	PHF19 (isoform b)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Preliminary experiments in Human Kidney, Skeletal Muscle and Placenta lysates gave no specific signal but low background (at antibody concentration up to 1µg/ml). We would appreciate any feedback from people in the field - have any results been reported w
Protein Sequence	DQRLPSQSLSSKQKQH
Clonality	Polyclonal
MW	22.5
Entrez	26147
Dilution Range	ELISA: 1:128000, WB: 1 µg/ml
Expiration Date	12 months from date of receipt.