


Product Datasheet

STXBP3 antibody (orb195334)

Description	Goat polyclonal antibody to STXBP3
Species/Host	Goat
Reactivity	Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	MUNC18C
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 70kDa band observed in Mouse and Rat Spleen lysates (calculated MW of 68.0kDa according to Rat NP_446089.1). Recommended concentration: 0.3-1µg/ml.
Protein Sequence	SKSIRRCKEINIS
Clonality	Polyclonal
MW	67.8
Entrez	6814
Dilution Range	ELISA: 1:128000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse (A) and R...