


Product Datasheet

STXBP3 antibody (orb195335)

Description	Goat polyclonal antibody to STXBP3
Species/Host	Goat
Reactivity	Human, Porcine
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	MUNC18C (aa514-527)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 60kDa band observed in Human Platelets lysates (calculated MW of 67.8kDa according to NP_009200.2). Recommended concentration: 0.1-0.3µg/ml.
Protein Sequence	RQKPRANYLED RKN
Clonality	Polyclonal
MW	67.8
Entrez	6814
Dilution Range	ELISA: 1:128000, WB: 0.1-0.3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Platelets...