


Product Datasheet

AMACR antibody (orb195338)

Description	Goat polyclonal antibody to AMACR
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	AMACR
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: pprox 40-42kDa band observed in Human Kidney and Liver lysates, and in lysates of cell line CaCo-2 and approx 40-45kDa band observed in lysates of cell line LNCaP (calculated MW of 42.4kDa according to NP_055139.4) and 43.9kDa according to NP_001161067.1
Protein Sequence	QLNSDKIIESN
Clonality	Polyclonal
MW	42.4; 22.2; 43.9
Entrez	23600
Dilution Range	ELISA: 1:64000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of CaCo-2 (A) and ...



Western blot analysis of Human Kidney (A) and ...