



Product Datasheet

AMACR antibody (orb195339)


Description	Goat polyclonal antibody to AMACR
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	AMACR (aa312-326)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 40-42kDa band observed in Human Kidney, Liver and Colorectal cancer lysates, and in lysates of cell line CaCo-2, and approx 40.45kDa band observed in lysates of cell line LNCaP (calculated MW of 42.4kDa according to NP_055139.4) and 43.9kDa accordi. IHC: Paraffin embedded Human Kidney. Recommended concentration: 5µg/ml.
Protein Sequence	RGSFITSEEQDVSPR
Clonality	Polyclonal
MW	42.4; 22.2; 43.9
Entrez	23600
Dilution Range	ELISA: 1:64000, WB: 0.1-0.3 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Human Ki...



Western blot analysis of CaCo-2 (A) and ...



Western blot analysis of Human Kidney (A...