

Product Datasheet

ITGB2 antibody (orb195343)

Description	Goat polyclonal antibody to ITGB2
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	ITGB2 / LFA-1 (aa86-99)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 90kDa band observed in lysates of cell lines K562 and U937 (calculated MW of 84.8kDa according to NP_00202.2). Recommended concentration: 0.3-1µg/ml. Primary incubation was 1 hour.. IHC: Paraffin embedded Human Spleen. Recommended concentration: 5µg/ml.
Protein Sequence	ETQEDHNGGQKQLS
Clonality	Polyclonal
MW	84.8
Entrez	3689
Dilution Range	ELISA: 1:128000, WB: 0.3-1 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of K562 (Lane 1) a...