


Product Datasheet

SP140 antibody (orb195347)

Description	Goat polyclonal antibody to SP140
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	SP140 (aa29-41)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 140kDa band observed in Human Spleen lysates (calculated MW of 98.2kDa according to NP_009168.4). Recommended concentration: 1-3µg/ml.
Protein Sequence	EGQNLQEQVCPEP
Clonality	Polyclonal
MW	98.2; 20.0; 95.3; 91.7; 86.0;
Entrez	11262
Dilution Range	ELISA: 1:128000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Spleen ly...