


Product Datasheet

GNL3 antibody (orb195355)

Description	Goat polyclonal antibody to GNL3
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	GNL3 (aa115-126)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	This antibody is expected to recognize both reported isoforms (NP_055181.3; NP_996561.1). Reported variants represent identical protein: NP_996561.1, NP_996562.1.WB: Approx 65kDa band observed in nuclear lysates of cell line K562 (calculated MW of 62.0kDa according to NP_055181.3). Recommended concentration: 0.03-0.1µg/ml.
Protein Sequence	TENKAKSGKQNS
Clonality	Polyclonal
MW	62.0; 60.5
Entrez	26354
Dilution Range	ELISA: 1:128000, WB: 0.03-0.1 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of K562 nuclear ly...



Immunohistochemical staining of Human Sm...