


Product Datasheet

CD74 antibody (orb195359)

Description	Goat polyclonal antibody to CD74
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	CD74 (aa159-171)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	This antibody is expected to recognize isoform a (NP_001020330.1) and isoform b (NP_004346.1).WB: Approx 33kDa band observed in Human Tonsil lysates (calculated MW of 33.5kDa according to NP_001020330.1). Recommended concentration: 0.01-0.03µg/ml.
Protein Sequence	KGSFPENLRHLKN
Clonality	Polyclonal
MW	33.5; 26.4
Entrez	972
Dilution Range	ELISA: 1:128000, WB: 0.01-0.03 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Tonsil ly...