


Product Datasheet

TNS1 antibody (orb195362)

Description	Goat polyclonal antibody to TNS1
Species/Host	Goat
Reactivity	Bovine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Target	Tensin 1 (aa902-914)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	IHC: In paraffin embedded Human Kidney shows preferential staining in cytoplasm of PCT over DCT. Recommended concentration, 2-4 µg/ml.
Protein Sequence	QYENQSPEATSPR
Clonality	Polyclonal
MW	186
Entrez	7145
Dilution Range	ELISA: 1:32000, IHC-P: 2-4µg/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Human Kidney