


Product Datasheet

NRXN1 antibody (orb19537)


Description	Goat polyclonal antibody to NRXN1
Species/Host	Goat
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, WB
Target	Neurexin 1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 150kDa band observed in Human Brain (Cerebellum) lysates (calculated MW of 162Da according to NP_004792.1). Recommended concentration: 0.05-0.1 µg/ml.
Protein Sequence	KEKQPSSAKSSNKN
Clonality	Polyclonal
MW	162; 46.9
Entrez	9378
Dilution Range	ELISA: 1:128000, WB: 0.5-1 µg/ml, IHC-P: 10ug/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Brain (Ce...)



Western blot analysis of HEK293 lysate U...



Flow cytometric analysis of paraformalde...