



Product Datasheet

Ogg1 antibody (orb19538)

Description	Goat polyclonal to Ogg1
Species/Host	Goat
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	Ogg1 (mouse)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:128000.WB: Approx 37kDa band observed in Rat Spleen lysates (calculated MW of 38.7kDa according to Rat NP_110497.1 and 38.8kDa according to Mouse NP_035087.2). Recommended concentration: 0.1-0.3 µg/ml.
Protein Sequence	HPKTSQAKGPSPLANK
Clonality	Polyclonal
MW	38.8, 38.8, 45.8, 47.2, 39.7, 22.2, 40.1, 36.1
Dilution Range	ELISA: 1:128000, WB: 0.1-0.3 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Rat Spleen lysa...



Immunohistochemical staining of Human Ki...