



Product Datasheet

PDE11A antibody (orb19541)

Description	Goat polyclonal antibody to PDE11A
Species/Host	Goat
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	PDE11A
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:2000.WB: Approx 60-65kDa band observed in Human Skeletal Muscle and Human Prostate lysates (calculated MW of 65.7kDa according to NP_001070826.1). Recommended concentration: 0.5-1.5 µg/ml.
Protein Sequence	NSFKESMEKSSYSD
Clonality	Polyclonal
MW	105, 78.1, 65.7, 55.7
Entrez	50940
Dilution Range	ELISA: 1:2000, WB: 0.5-1.5 µg/ml, IHC-P: 10 µg/ml
Expiration Date	12 months from date of receipt.



orb19541
(1µg/ml)
staining of
Human Pros...



orb19541
(10µg/ml)
staining of
paraffin...