



Product Datasheet

XRCC6 antibody (orb19544)

Description	Goat polyclonal antibody to XRCC6
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	KU70 / XRCC6
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:2000.WB: Approx 70kDa band observed in lysates of cell lines A431, HeLa and HepG2 (calculated MW of 69.8kDa according to NP_001460.1). In transfected HEK293 transiently expressing KU70 a band of approx. 70kDa is observed. This band is not observed in the non-transfected HEK293.Recommended concentration: 0.1-0.3 µg/ml.
Protein Sequence	YKTEGDDEEAAAAEQ
Clonality	Polyclonal
MW	69.8
Entrez	2547
Dilution Range	ELISA: 1:2000, WB: 0.1-0.3 µg/ml
Expiration Date	12 months from date of receipt.



orb19544
(0.1ug/ml)
staining of
HeLa cel...



HEK293
overexpressing
KU70
(RC204048)
an...