



Product Datasheet

APH1A antibody (orb19545)


Description	Goat polyclonal antibody to APH1A
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	APH1A
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:128000.WB: Approx 26kDa band observed in Human Brain (Cerebral Cortex) lysates (calculated MW of 26.8kDa according to NP_057106.2). Recommended concentration: 0.05-0.1 µg/ml.
Protein Sequence	HVTDRSDARLQYG
Clonality	Polyclonal
MW	28.9, 26.8
Entrez	51107
Dilution Range	ELISA: 1:128000, WB: 0.05-0.1 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



orb19545 (0.05ug/ml)
staining of Human
B...



Immunohistochemical
staining of Human
Ki...



Immunohistochemical
staining of Human
Pl...