


Product Datasheet

IL12 p40 antibody (orb19546)


Description	Goat polyclonal antibody to IL12 p40
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	IL12B / IL12p40
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:16000.IHC: In paraffin embedded Human Colon shows cytoplasm staining of individual cells between columns. In paraffin embedded Human Kidney shows cytoplasm staining of epithelial cells of some DCT. Recommended concentration, 4-6ug/ml.WB: Approx 37kDa band observed in Human Liver lysates (calculated MW of 37.2kDa according to NP_002178.2). Recommended concentration: 0.1-0.3 µg/ml.
Protein Sequence	QGKSKREKKDRVFTD
Clonality	Polyclonal
MW	37.2
Entrez	3593
Dilution Range	ELISA: 1:128000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



orb19546
(0.1ug/ml)
staining of
Human Li...



orb19546
(4ug/ml)
staining of
paraffin e...



Western blot
analysis of
Mouse Liver
(la...