



Product Datasheet

MYB antibody (orb19550)


Description	Goat polyclonal antibody to MYB
Species/Host	Goat
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	c-Myb (aa281-294)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:8000.IHC: In paraffin embedded Human Kidney shows cytoplasm staining of renal tubules and in Human Breast cytoplasm staining of ductal epithelial cells. Recommended concentration, 3-6ug/ml.WB: Approx 75kDa band observed in Human Spleen lysates (calculated MW of 72.3kDa according to NP_005366.2). Recommended concentration: 0.3-1 µg/ml.
Protein Sequence	QRHYNDEDPEKEKR
Clonality	Polyclonal
MW	72.3
Entrez	4602
Dilution Range	ELISA: 1:4000, WB: 1-3 µg/ml, IHC-P: 4 µg/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Human Ki...



Immunohistochemical staining of Human Ki...



Western blot analysis of MOLT4 cell lysa...