

Product Datasheet

RAB5A antibody (orb19552)

Description

RAB5A antibody

Species/Host

Rabbit

Reactivity

Guinea pig, Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ICC, IF, IHC-P, WB

Applications
Immunogen

KLH conjugated synthetic peptide derived from human RAB5A. Please contact us for the exact immunogen sequence. The peptide is available as orb393975.

Target

RAB5A

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW



24 kDa

Uniprot ID
P10683
NCBI
015973, 057057
Entrez
51083
Dilution Range


IF/ICC: 1:100-500, IHC-P: 1:100-500, WB: 1:200-2000

Expiration Date


12 months from date of receipt.



Immunohistochemical staining of paraffin...



IF analysis of guinea pig spleen tissue ...



IHC-P staining of guinea pig spleen tiss...