


Product Datasheet

TAIP-12 antibody (orb19553)

Description	Goat polyclonal to TAIP-12
Species/Host	Goat
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	TAIP-12 / FAM130A1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:64000.WB: Approx 60kDa band observed in Human Heart and Mouse Heart lysates (calculated MW of 58.5kDa according to NP_700456.1). Recommended concentration: 0.1-0.3 µg/ml.
Protein Sequence	KRQKQLRRKNVRFDQ
Clonality	Polyclonal
MW	58.5
Entrez	81566
Dilution Range	ELISA: 1:64000, WB: 0.1-0.3 µg/ml
Expiration Date	12 months from date of receipt.



orb19553
(0.1ug/ml)
staining of
Human He...