

Product Datasheet

PAX6 antibody (orb19559)

Description	Goat polyclonal antibody to PAX6
Species/Host	Goat
Reactivity	Bovine, Frog, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, WB
Target	PAX6 (internal)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 50kDa band observed in Mouse Eye lysates (calculated MW of 48.2kDa according to Human NP_001595.2, and to Mouse NP_038655.1). Recommended concentration: 0.3-1 µg/ml.
Protein Sequence	REEKLRNQRRQASN
Clonality	Polyclonal
MW	46.7; 48.2
Entrez	5080
Dilution Range	ELISA: 1:32000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



orb19559
(0.5ug/ml)
staining of
Mouse Ey...