

Product Datasheet

GAB2 antibody (orb19563)

Description	Goat polyclonal antibody to GAB2
Species/Host	Goat
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	GAB2
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:64000.WB: Approx 80kDa band observed in lysates of cell line K562 (calculated MW of 74.5kDa according to NP_536739.1). Recommended concentration: 0.5-1.5 µg/ml.
Protein Sequence	QVDKEKTQALQN
Clonality	Polyclonal
MW	74.5; 70.0
Entrez	9846
Dilution Range	ELISA: 1:64000, WB: 0.5-1.5 µg/ml
Expiration Date	12 months from date of receipt.



orb19563
(0.5ug/ml)
staining of
K562 cel...