



Product Datasheet

TPD52L2 antibody (orb19569)


Description	Goat polyclonal antibody to TPD52L2
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	TPD52L2 / D54
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:16000.WB: Approx 28kDa band observed in Human Breast Cancer and Human Brain lysates (calculated MW of 24.9kDa according to NP_955392.1). Recommended concentration: 0.1-0.3 µg/ml. See also our D52 and D53 products.
Protein Sequence	SGDKPLSDPAP
Clonality	Polyclonal
MW	24.9; 22.5; 23.8; 21.4; 22.2; 19.9
Entrez	7165
Dilution Range	ELISA: 1:16000, WB: 0.1-0.3 µg/ml, IHC-P: 3.75 µg/ml
Expiration Date	12 months from date of receipt.



orb19569 (0.1ug/ml)
staining of Human
Br...



Immunohistochemical
staining of Human
Co...



Immunohistochemical
staining of Human
To...