


Product Datasheet

THRA antibody (orb19580)


Description	Goat polyclonal antibody to THRA
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	THRA
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 60+50kDa bands observed in Human Bone Marrow lysates (calculated MW of 54.8kDa according to NP_003241.2 and of 46.8kDa according to NP_955366.1). Recommended concentration: 0.5-1,5 µg/ml
Protein Sequence	KEEMIRSLQQRP
Clonality	Polyclonal
MW	46.8; 54.8
Entrez	7067
Dilution Range	ELISA: 1:32000, WB: 0.5-1ug/ml, IHC-P: 3.75 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Bone Marrow lysates



Immunohistochemical staining of Human Bone Marrow tissue



Immunohistochemical staining of Human Brain tissue