

Product Datasheet

NEDD9 antibody (orb19581)

Description	Goat polyclonal to NEDD9
Species/Host	Goat
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	NEDD9
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	This antibody is expected to recognise isoform 1 (NP_006394.1) only.ELISA: Peptide ELISA: antibody detection limit dilution 1:8000.WB: Approx 100kDa band observed in Human Kidney, Lung and Placenta lysates (calculated MW of 92.9kDa according to NP_006394.1). Recommended concentration: 0.1-0.3 µg/ml.
Protein Sequence	NKPQNKDDLDR
Clonality	Polyclonal
MW	92.9
Entrez	4739
Dilution Range	ELISA: 1:8000, WB: 0.1-0.3 µg/ml
Expiration Date	12 months from date of receipt.



orb19581
(0.1ug/ml)
staining of
Human Ki...