

Product Datasheet

C4BPA antibody (orb19586)

Description

C4BPA antibody

Species/Host

Rabbit

Reactivity

Guinea pig, Human, Mouse, Rat

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications

Immunogen

KLH conjugated synthetic peptide derived from human C4BPA. Please contact us for the exact immunogen sequence. The peptide is available as orb393977.

Target

C4BPA

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW

67 kDa

Uniprot ID

[P10683](#)

NCBI

[015973](#), [057057](#)

Entrez



[51083](#)

Dilution Range


IHC-P: 1:100-800, WB: 1:200-1000

Expiration Date


12 months from date of receipt.



Immunohistochemical staining of guinea p...



WB analysis of guinea pig lymph (lane 1)...



rat stomach (lane 4)