



Product Datasheet

EGR2 antibody (orb19587)


Description	Goat polyclonal antibody to EGR2
Species/Host	Goat
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	EGR2
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:64000.IHC: In paraffin embedded Human Colon shows nuclear staining of absorptive cells. Recommended concentration, 3-5ug/ml.WB: Approx 55kDa band observed in lysates of cell line HepG2 (calculated MW of 50.3kDa according to NP_000390.2). Recommended concentration: 0.03-0.1 µg/ml.
Protein Sequence	HGTAGPDRKPFPC
Clonality	Polyclonal
MW	50.3
Entrez	1959
Dilution Range	ELISA: 1:8000, WB: 2-6 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of HepG2 lysate us...



Immunohistochemical staining of Human Co...



Immunohistochemical staining of Human Sm...